

# myGaze® Assistive 2

GAZE CONTROL MADE EASY AND AFFORDABLE



[WWW.MYGAZE.COM](http://WWW.MYGAZE.COM)

# Designed for Simplicity and Ease of Use

## EYEMOUSE PLAY

World's first  
easy-to-use Eye Mouse



## STATE-OF-THE-ART MYGAZE EYE TRACKER

Tracks where your eyes are looking  
and turns that into computer control

## PLUG & PLAY

No more expensive training  
and support!

## NO-FAIL START

0- and 1-point calibration for  
early users

## UNIQUE CLASSROOM FEATURES

Designed with students and  
teachers in mind

## AS YOUR SKILLS IMPROVE ...


you find more possibilities  
within the same myGaze.  
No need to upgrade!

# Designed for Simplicity and Ease of Use

## WHO IS IT FOR?

### CLASSROOM

Designed to help special need students at any ability learn, participate and communicate



### AAC & PC CONTROL

Provides pixel accuracy and independence for advanced gaze users

## WHAT CAN YOU DO WITH IT?

Everything you can do with a mouse! The myGaze Assistive System will follow your gaze on the screen and let you move the cursor and issue mouse commands through multiple input methods such as “dwell”, “blink” and “switch”.

**THE ONLY EYE GAZE  
SYSTEM FOR EVERY STEP  
IN A CHILD'S PROGRESS**

**AN EASY START**  
No more fear of failure or discouragement! With our *calibration-free mode* every child can get started at their own ability level and pace.


**INCLUSION  
IN CLASSROOM**  
With our unique set of features and bundles teachers can *asses*, *include* and *engage* their students to achieve high educational goals.

**BUILDING  
INDEPENDENCE**  
As you improve your skills myGaze allows you to do ever more: choose your own settings, open desktop applications, participate and communicate with peers.

WWW.MYGAZE.COM

# The Eye Gaze for the Classroom

Award winning gaze control for every level of ability



## PLUG & PLAY

Easy 3-step set up, no-fail calibration and you are ready to go! No expensive training & support required.



## GAZE ANALYSIS REPORT

Heat map, gaze trace line, and playback video support the teacher's understanding with every activity in the myGaze bundles of Inclusive Technology.




## MYGAZE WITH SWITCH

With our EyeMouse, Switch users can now enjoy so many more 'cursor' activities on the PC.



## NON-INTRUSIVE TEACHER CONTROL

Keyboard shortcuts, mouse-over-gaze control ... A full set of classroom features makes the teacher-student interaction easy and natural.



## INTUITIVE POSITIONING GUIDE

The 'glasses' positioning is fun and easy even for young users. The color cues provide live feedback at all time!



## FEATURES FOR ADVANCED USERS

Look at camera to call the menu, 9-point calibration for pixel accuracy and many customization options in our fully gaze-accessible menu.

## CLASSROOM BUNDLES





## For communication, PC control and more with myGaze and Grid 3

The new Grid Pad Eye with myGaze system provides an innovation solution which is simple to use, empowering you to communicate, control your PC and your environment.



Our reliable and accurate gaze control, combined with the powerful Grid 3 software from Smartbox, is suitable for everyone. Whether you are just starting out with eye gaze, and beginning to learn early access skills, or an experienced eye gaze user, looking for fast and independent communication, there are tools and resources in the myGaze + Grid 3 solution to suit your needs.



Symbol  
Communication



Text  
Communication



Interactive  
Learning



Accessible  
Apps



Environment  
Control



Computer  
Control

## Stories: Hear the Users

**Declan** is an 8-year-old kid, who like many other has *substituted his Switch for a myGaze*. By translating his eyes into a cursor on the screen, myGaze gives him the opportunity to study with the same resources as his peers, participate more in class and develop his communication skills.

*“With the unique classroom features and full teacher control, myGaze has made my interaction with Declan very natural and non-intrusive.”*

By **Jennette Holden**, Pendle View Primary School, therapist of Declan

**Sally** has a *degenerative condition*, which had made it difficult to involve her in class. She very quickly understood she was activating the program with her eyes and this was apparent through her eye movements on the screen.

*“This equipment is a vital tool for educators. It allows a greater insight into their awareness and responses.”*

By **Carol Allen**, ICT and SEN advisor in North Tyneside, therapist of Sally

**Sonja**, the girl with *Rett Syndrome*, is able to make choices, communicate her needs and operate simple PC applications, all within a 3-month trial with myGaze.

*“With myGaze device we see our daughter make choices that she has never had the ability to make.”*

By **Sara Barbee**, mother of Sonja, Portland, US

DESIGNED FOR SIMPLICITY AND EASE OF USE

[INFO@MYGAZE.COM](mailto:INFO@MYGAZE.COM)

Or ask your preferred reseller or therapist, worldwide.

[WWW.MYGAZE.COM](http://WWW.MYGAZE.COM)

